

WE MOVE YOU



Andago®

QUICK REFERENCE GUIDE



Make sure you have read the user manual and watched the instructional video before using the Andago V2.0 with patients. The video is available on knowledge.hocoma.com


PREPARATION

1. Place harness on patient (check symmetry)
2. Turn Andago on (next to rear right bumper)
3. Write down patient ID (from screen)



PATIENT SETUP

1. Lift patient (until standing)
2. Adjust body weight support
3. Select default fall limit (hold button)



Holding the override button activates joystick at any time

OPTIONAL STEPS Patient Setup

- 1a. Select asymmetrical body weight support
- 3a. Adjust fall limit using patient lift buttons

TRAINING

1. Select training mode (patient following / straight line / manual)
2. Select speed limit
3. Start training

Training

- 3a. Adjust turning offset

FINISH SESSION

1. Stop training
2. Lower patient
3. Press Reset button 2 sec (to save data)
4. Open harness and bring patient out of Andago

